ACT
Adjustable Continence Therapy for Women
Uromedica, Inc.

Instructions for Use

Please read all information carefully.

All first-time implanting surgeons must have received an in-service training session with a surgeon previously qualified on the ACT implantation procedure.

Failure to properly follow instructions may result in improper functioning of the device and lead to injury.

Important: This package insert is designed to provide instructions for use of the Uromedica ACT single use device and Uromedica reusable surgical instruments. **It is not a comprehensive reference to surgical technique for correcting Stress Urinary Incontinence (SUI).** The device should be used only by physicians trained in the surgical treatment of SUI. These instructions are recommended for general use of the device. Variations in use may occur in specific procedures due to individual technique and patient anatomy.

I. Description of ACT System

The ACT therapy consists of:

- A periurethral implantable prosthesis which has a two-lumen extruded tube (2.4 mm outside diameter) with an expandable balloon on the end that is implanted adjacent to the bladder neck. The balloon contains a 1mm radiopaque marker positioned 1mm from the distal tip. Balloon diameter is 25 mm at maximum fill volume (8 mL). The implant contains a subcutaneous titanium port on the end of the tube that is used by the physician for post-operative adjustment.
- A pushwire of 0.8 mm OD and 15 cm in length with a flexible tip on which the ACT device is mounted and is used to deliver the ACT through an Implantation Instrument Set.
- A syringe with a non-coring needle.
- Material exposure to user (Implant): 100% silicone.
 - The maximum surface area exposure of this product includes two 9 cm devices, each at maximum fill volume (8 mL). This maximum surface area exposure is 7,802 in².
 - The maximum mass of this product includes two 9 cm devices, each at maximum fill volume (8 mL). This maximum mass is 17.8 g.
- Maximum fill volume of isotonic contrast solution per device: 8 mL.

II. Contents of ACT Patient Package:

Table 1 – ACT Model Numbers

Model Number	Device Length	Description
800016-01	5 cm	ACT Patient Kit: Two (2) prosthetic implants
800016-02	6 cm	ACT Patient Kit: Two (2) prosthetic implants
800016-03	7 cm	ACT Patient Kit: Two (2) prosthetic implants
800016-04	8 cm	ACT Patient Kit: Two (2) prosthetic implants
800016-05	9 cm	ACT Patient Kit: Two (2) prosthetic implants
800021-01	5 cm	ACT Revision Kit: One (1) prosthetic implant
800021-02	6 cm	ACT Revision Kit: One (1) prosthetic implant
800021-03	7 cm	ACT Revision Kit: One (1) prosthetic implant
800021-04	8 cm	ACT Revision Kit: One (1) prosthetic implant
800021-05	9 cm	ACT Revision Kit: One (1) prosthetic implant

All configurations also include:

One (1) 0.60 mm (23G) x 28.6 mm Huber Needle

One (1) 5 mL Syringe with Luer Lock Connector
One (1) Spare Pushwire
One (1) Instructions for Use Insert
One (1) Product Label Set

III. Uromedica Surgical Instruments (sold separately)

- A non-sterile Implantation Instrument Set (P/N 750009) consisting of a U-channel sheath and sharp-tipped trocar.
- A non-sterile Blunt Trocar (P/N 750022)
- A non-sterile Tissue Expanding Device (TED II) (P/N 750034)

The Uromedica Blunt Trocar and Uromedica TED are both intended to be used with the U-channel sheath of the Uromedica Implantation Instrument Set. The U-channel sheaths mates with these devices during use, per the instructions below.

IV. Intended Purpose and Indication

Intended Purpose

The ACT prosthetic implant device is intended, via initial and post-operative volume adjustments, to increase urethral coaptation at the bladder neck to eliminate unwanted urine loss.

Indication

The ACT system is a permanent implant designed for the correction of SUI in adult female patients.

V. Contraindications

- Patients with active systemic or urinary tract infections.
- Patients who are pregnant.
- Patients with unmanageable detrusor instability.
- Patients with reduced bladder compliance.
- Patients with pretreatment urge incontinence symptoms refractory to medication.
- Patients with neurogenic bladder conditions including detrusor sphincter dyssynergia (DSD)
- Patients with residual volume greater than 100 mL after voiding.
- Patients who have received, are presently receiving or will receive pelvic-region radiotherapy within 6 months of ACT system implantation.
- Patients with bladder cancer or bladder stones.
- Patients with bleeding disorders.

VI. Warnings

- The device is an implantable device and is intended for a single use only. Do not re-sterilize. Re-sterilization could result in the mechanical failure of the device and place the patient at undue risk of procedural failure. Placement of the ACT device must be performed by a physician.
- If the physician chooses to use X-ray, adequate attention must be given to the potential radiation exposure.
- Perforation of the bladder, urethra, or vagina during dilation of the tissue from the labia to the bladder neck should be avoided.
- Do not use excessive force to advance or withdraw the device if resistance is encountered.
- Do not attempt to place or introduce the device without the supplied pushwire.
- Do not place the device in a location where the patient may have problems with urinary retention or erosion.
- Do not inflate the balloon to more than 8.0 mL.
- Use the correct mixture of sterile water and contrast media to provide an isotonic contrast solution for inflation of the ACT balloons. See attached reference chart (Appendix A).
- Do not use any of the non-sterile accessory surgical instruments without prior cleaning and steam sterilization as indicated in section IX.
- Do not place the device in the urethral mucosa. This could result in balloon erosion or migration.
- Device wear or puncture may result if the device comes in contact with itself or another device, or from surgical staples, steel suture, or other artifacts from previous interventions present in the region where the device will

- be positioned. The physician is cautioned to carefully evaluate the implant location.
- The ACT system is MR-Conditional. A patient with this device can be scanned safely in an MR (Magnetic Resonance Imaging) system immediately after placement under certain conditions as detailed in the MRI Information section provided in these Instructions for Use (IFU).
- Patients who have had previous surgeries for incontinence may have a greater incidence of complications than patients who have not had previous surgeries for incontinence.
- Device leak may result due to long-term calcification on the surface of the balloon.
- Resolution of adverse events (e.g. device erosion, prosthetic infection, device wear) or suspected device malfunction or change in performance may require explant of one or both of the devices. The explanted device(s) can be replaced in 60 to 90 days.
- Patients who are sexually active may have a greater incidence of complications than patients who are not.

VII. Precautions

- Do not use the ACT system prior to reading and understanding the Instructions for Use. Product packaging that is damaged or unintentionally opened shall not be used.
- Do not attempt to re-sterilize the ACT device.
- To reduce the risk of needle sticks, always keep the plastic needle cap on the needle until the needle is mated with the syringe and only remove the needle cap immediately prior to using the needle.
 - Proper Disposal: do not attempt to recap a used needle. Do not attempt to bend or break the needle. Immediately dispose of the entire unit (needle, syringe body, and cap) into a puncture-resistant biohazard sharps container.
 - Do not reach into sharps container.
 - Do not leave used needle on the table, even if the cap is fully covering needle.
- Explanted devices or devices potentially contaminated should be handled as biohazardous waste.
- The ACT system should only be used by a physician trained 1) in either urologic or urogynecologic surgical procedures and 2) who has also completed the in-service training session on the ACT implantation technique conducted by a previously qualified surgeon.
- Patients undergoing a surgical implantation of the ACT device should be pre-medicated with an antibiotic to reduce the possibility of a postoperative infection.
- Once placed, the ACT device should be viewed with vaginal ultrasound, fluoroscopy, or digitally palpated from the anteriolateral vaginal wall.
- It is the physicians' responsibility to determine the correct length of the ACT device.
- Do not implant the ACT device without use of the Uromedica, Inc. Implantation Instrument Set and the provided non-coring needle for filling the balloons.
- Do not implant just one ACT device EXCEPT in the case of bladder, urethra or vaginal wall perforation. If a perforation occurs, do not implant an ACT device on the affected side for at least 60 days. An effective procedure calls for two devices to be placed.
- Improper placement of the ACT device may lead to a sub-optimal patient clinical response or patient discomfort.
- Maintain the sterile field throughout the procedure.
- If the patient is to undergo any further examination or surgery in the pelvic region after the ACT device has been implanted, the physician must consider possible adverse effects.
- Do not handle the ACT device, particularly the balloon, with any sharp instruments or lay it against an abrasive material. Failure of the device may result. Do not grasp the port or strain relief with clamps.
- If the patient is allergic to contrast media, the ACT balloons should be filled with sterile saline instead of isotonic contrast solution.
- Use the same type of contrast media as used during the original filling of the balloons when performing a post surgical adjustment.

- Postoperatively, the patient is recommended to refrain from heavy lifting and/or exercise for 3 to 4 weeks and intercourse for one month.
- Should bleeding or other problems occur, the patient is advised to contact the surgeon immediately.
- Patients with a current clinical indication of cystocele, enterocele, rectocele urethrocele or uterine prolapse may have a lower rate of efficacy and should have these repaired prior to implant of the ACT system.
- The Tissue Expanding Device II should not be rotated while the tip is extended. Rotating the tool while the tip is extended may increase the risk of unintended tissue/organ damage.
- The tip of the Tissue Expanding Device II should not be extended towards of the urethra. Extending the tip towards the urethra may increase the risk of perforation.
- The Tissue Expanding Device II should not be used to grasp material. Grasping material may damage the tool.
- Patients with intrinsic sphincter deficiency (ISD) from neurogenic causes may have a worse functional result and may have a higher risk for intra-operative and post-operative complications.

If the patient becomes pregnant subsequent to implantation of the ACT device, the primary physician must closely monitor the patient for any adverse effects associated with the device during the pregnancy.

VIII. Potential Adverse Events (Residual Risks)

Potential adverse events from the ACT system, the attendant equipment and supplies, or the implant procedure include:

Allergic Response
 Anesthetic risks
 Device Calcification
 Device Malfunction
 Device Wear
 Erosion of the device at the bladder, labia, urethra, or vagina
 Explant of one or both devices
 False Channel Creation
 Hydronephrosis
 Intraoperative perforation (e.g. bladder, urethra, vagina)
 Pain or Discomfort
 Procedure/adjustment not able to be completed
 Prosthetic Infection
 Prosthetic Migration
 Urinary Complications (e.g. Urinary Difficulty, Frequency, or Urgency, Worsening Incontinence*)
 Urinary Retention
 Urinary Tract Infection
 Vascular Complications (e.g. Bleeding, Hematoma, Induration, Inflammation)
 Wound Dehiscence
 Wound Infection

* "Worsening Incontinence" includes subjective patient-perceived worsening.

**Any serious incident associated with the device must be promptly reported to the manufacturer and the relevant competent authority in the Member State where the user or patient is located.

IX. Sterility and Storage

The contents of the undamaged and sealed packaged ACT device are STERILE. The package has been sterilized with gamma irradiation.

Store product at room temperature with a temperature limit of -30°C to 60°C. Keep the product dry.

The Uromedica Implantation Instrument Set, Blunt Trocar, and Tissue Expanding Device II (TED II) are supplied non-sterile and must be sterilized before initial use and cleaned and sterilized after each use. Follow the reprocessing instructions provided in the Instructions for Use provided with the surgical instruments.

X. MRI Information

MRI Safety Information	The ACT, Adjustable Continence Therapy for Women is MR Conditional. A patient with the ACT,
-------------------------------	---

	Adjustable Continence Therapy for Women may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.
Parameter	Condition of Use/Information
Static Magnetic Field Strength (T)	1.5-Tesla or 3.0-Tesla
Static Magnetic Field Orientation	Horizontal
Maximum Spatial Field Gradient (T/m and gauss/cm)	50-T/m (5,000-gauss/cm)
RF Excitation Polarization	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil	Any transmit RF coil may be used
Receive-Only RF Coil	Any receive-only RF coil may be used
MR System Operating Mode	Normal Operating Mode
Maximum Whole Body Averaged SAR (W/kg)	2-W/kg (Normal Operation Mode)
Scan Duration and Wait Time	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks)
MR Image Artifact	The presence of this implant produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant is located in the area of interest.

XI. Instructions for Implantation

Equipment Requirements

- Image Intensifier and sterile cover (optional)
- Camera and light source

Instrument Requirements

- Rigid or flexible cystoscopy set up
- Regular Surgical instrument set
- No 10 blade and handle
- Debaquey Forceps
- Rubber shod plain forceps
- Adson Forceps
- Cairns Tissue Forceps
- McIndoe / Metzenbaum Scissors
- Dressing Scissors
- 7" Needleholder
- Uromedica Implantation Instrument Set (sharp trocar and U-channel sheath)
- Uromedica Blunt Trocar
- Uromedica Tissue Expanding Device (TED II)

Accessories

- Jug
- Gallipots x2
- Kidney dish x2
- Bladder syringe
- Syringe
- 3/0 Vicryl rb needle
- 4/0 skin suture
- Foley catheter (typically 18Fr) and drainage bag

Solutions

- One- liter normal sterile saline (0.9%) or sterile water for irrigation during cystoscopy.
- Isotonic contrast solution for filling the device. See Appendix A.
- Sterile saline for soaking the ACT device prior to implantation.
- Use 50 - 100 mls of contrast and dilute with normal sterile saline (0.9%) to reach a total volume of 250 mls.

Preparing the Patient

- Start oral antibiotic therapy.
- Sedatives may be given preoperatively.
- Place the patient in the lithotomy position.
- Shave and prepare the perineum, proximal/anterior and medial thigh and vulva with an iodophor (if the patient is not allergic) or equivalent surgical preparation.
- Perform formal vaginal preparation
- Place an appropriately sized (typically 18 French) Foley catheter. Fill the bladder with 200 ml of dilute contrast. Fill the Foley catheter balloon with 10 ml of sterile water or with contrast to better visualize the bladder neck when intraoperative x-ray is used.

Administering Anesthetic

- Administer either general, spinal, or local anesthetic which should be the preference of the surgeon.
- When administering local anesthetic, deliver along the entire tissue path through which the Implantation Instrument Set will be passed.

Preparing the Contrast Media

- Prepare an isotonic filling solution by diluting a contrast media chosen from Appendix "A" with the prescribed amount of sterile water. The full volumes indicated need not be mixed but the ratios should remain the same so that water does not pass into or out of the device over time due to osmosis. Larger volumes facilitate more accurate mixing.

Prepare the Uromedica Re-usable Surgical Instruments

- The Uromedica Implantation Instrument Set, Blunt Trocar, and Tissue Expanding Device II (TED II) are supplied non-sterile and must be sterilized before initial use and cleaned and sterilized after each use. Follow the reprocessing instructions provided in the Instructions for Use provided with the surgical instruments.

ACT System Implantation

- Make a small incision between the labium majus and minus.
- The Instrument Implantation Set is guided by twisting it in a rotating motion under digital vaginal control in a lateroposterior direction to the pubic bone, just below the level of the urethral meatus. If there is scar tissue present, physicians often use the sharp tipped trocar that provides a cutting effect when twisted, facilitating guidance of the Implantation Instrument Set through fibrous tissue. The blunt tipped trocar can be used in place of the sharp tipped trocar, if preferred, for penetration of non-fibrous tissue. In order to ease the guidance of the Instrument Implantation Set, the designated route may be prepared using scissors.
NOTE: The Foley balloon should be pulled back into the bladder neck to allow for palpation of the Foley balloon providing a clear landmark of the bladder neck and implant location.
- In order to prevent unintended perforation of the vaginal wall during the procedure, make sure that the Instrument Implantation Set is run lateral to the full thickness of the vaginal wall, especially in patients with previous vaginal surgery. When the tip of the trocar reaches the pubic bone, the Instrument Implantation Set is then guided under digital control in a medioanterior direction. In order to reach the intrapelvic portion of the urethra, the pelvic diaphragm (M. transversus perinei profundus) must be penetrated. After the penetration of the pelvic diaphragm, the sharp trocar should be exchanged with the blunt tipped trocar to reduce risk of bladder or urethral perforation.
- When the tip of the trocar can be palpated at the urethrovaginal junction just below the bladder neck, remove the trocar. At this point, if the patient has fibrotic tissue, the Tissue Expanding Device (TED) may be used to dilate the tissue, creating a space for the ACT balloon. The TED can be locked within the U-Channel sheath. The TED handles are then squeezed together to open

the jaw, facilitating tissue dilation. The TED handles must be opened fully (which closes the jaw fully) prior to removing the TED from the U-Channel sheath.

5. Choose an ACT device with the appropriate length. An 8 cm length ACT device is chosen most often. In more obese women, a 9 cm device will often be needed. Centimeter markings on the trocar may aid in choosing the proper length ACT device. Once the tip of the trocar is at the incision site, count the number of centimeters marked on the back side of the trocar and add two centimeters for the final approximate length of the ACT device.
6. Place the needle (included) in the septum of the titanium port and fill the balloon with 1.0 to 2.0 ml of the previously prepared isotonic contrast solution. Aspirate two or three times to remove most of the air from the ACT device. Remove the needle.
7. Soak the ACT system in sterile saline (Solution C) until it is ready to be implanted.
8. Just before implantation, inflate the balloon with approximately 3.0 ml of the isotonic contrast solution to check for any leaks or excessive distortion of the balloon. Aspirate all the fluid from the balloon. It should form a tri-wing configuration to facilitate insertion. It may be necessary to form the balloon into a tri-wing configuration with your fingers as the fluid is aspirated. The balloon is then coated with surgical lubricant and one 'wing' is placed directly into the channel of the U-channel sheath.
9. Advance the ACT device with the balloon deflated through the U-channel sheath to the selected periurethral position using the pushwire. The position of the tip of the device can be visualized either by x-ray or transvaginal ultrasound and should be placed underneath the bladder neck.
10. Pull the U-channel sheath of the Implantation Instrument Set back about 2 cm. Penetrate the port septum with the needle (included) on the syringe and inflate the ACT balloon with 1.0 ml of the prepared isotonic contrast solution. Under x-ray visualization (or transvaginal ultrasound), the balloon should form a spherical shape. While making sure the balloon remains in the correct position, withdraw the U-channel sheath completely.
11. Check the position of the inflated ACT balloon digitally. There should be some mobility of the anterior vaginal wall against the balloon. This ensures that the balloon has not been inflated within tissue layers adjacent to the vaginal mucosa.
12. Leave the pushwire in place. Repeat the previous steps on the contralateral side unless performing a unilateral revision procedure.
13. Confirm symmetrical positioning of the ACT balloons with respect to the urethra either radiologically or by transvaginal ultrasound. Remove the pushwires from both ACT devices after balloon inflation.
14. Confirm the absence of urethral or bladder injury by cystoscopic examination.
15. Create a pocket in the subcutaneous tissue of the labium majus with scissors through the small skin incision in an anterior direction. Place the ACT port in this pocket. This procedure facilitates palpation of the port for adjustments of the ACT device in the future.
16. Proceed to create a pocket in the contralateral side in the same manner unless performing a unilateral revision procedure. Place the remaining ACT port in this pocket. Check that the tubing of the ACT devices is not kinked.
17. Close both incisions in two layers using resorbable suture for the subcutaneous tissue and the skin.
18. The patient should continue to take the same oral antibiotic given preoperatively for a maximum of 3 days.
19. Schedule the patient for a 6-week follow-up visit for assessment and adjustment as required.
20. Confirm the patient can void spontaneously. If the patient cannot void spontaneously, place a Foley catheter until anaesthesia effects have subsided, then have patient void. If the patient is still unable to void, decrease the balloon volume in 0.5 ml increments from both balloons until the patient can spontaneously void.

XII. Post-Operative Care

- The patient is usually discharged on the day of the surgery. Ensure that the patient can void spontaneously prior to discharge from the surgical facility. Most surgeons choose to send their patients home with a course of antibiotics. The patient should be advised on the use of absorbent pads until the desired level of continence is achieved. Adjustments may begin six weeks after the ACT implantation procedure. Inform the patient that multiple adjustments may be required to achieve the desired level of continence.
- Provide patient education about the operation. The patient should be advised to avoid undue strenuous activity immediately post-operatively as directed by the physician. Failure to do so could lead to device migration. Instruct the patient to contact the physician if they experience redness, swelling, or drainage of the incision site, experience a return of incontinence symptoms, or are unable to urinate. Provide the patient with information on the MRI conditional status of the ACT implant and the impact that ACT might have on travel security screenings (titanium and radiopaque components).
- After the postoperative healing period, the surgeon should maintain contact with the patient to evaluate the function of the device. During an evaluation, the physician should ask the patient about how ACT is functioning and if the patient has noticed any changes in the function. If the patient's incontinence symptoms have returned, or there is infection or erosion present, revision or removal surgery may be necessary. In the event of a revision surgery, follow the same explanation and implantation techniques outlined elsewhere in this document.

XIII. Postoperative Balloon Adjustment

- After edema from the surgery has subsided (approximately six weeks postoperatively) the patient may become incontinent, in which case, adjustments should be performed.
- Adjustments may be conducted in the office. Topical anesthetic may be used, if needed to minimize patient discomfort during palpation of the port and insertion of the needle.
- Do not increase the balloon volume until six weeks post-implantation. Increases in balloon volume performed before six weeks post-implantation are not advised and can result in device migration.
- Use a 0.60 mm (23G) x 28.6 mm needle for post operative balloon volume adjustments. Use of other needles may result in damage to the device port.
- Clean and prepare skin over the adjustment port. Palpate the implanted port to allow percutaneous adjustment of balloon volume by accessing the port using a 0.60 mm (23G) x 28.6 mm needle. The ACT balloon volume should be increased in 0.5 ml increments. To maintain isotonicity, use the same inflation media (diluted contrast media or sterile saline) when performing a balloon volume adjustment as was used during the initial inflation. CAUTION: If resistance is felt, do not continue pushing the needle into the port. Pull the needle back and ensure device tubing is not kinked. Pushing the needle into the edge of the port can cause damage to the port, which may require explant of the device.
- Make sure that patient can spontaneously void prior to leaving the office.
- If the patient experiences urinary retention, remove volume from each of the ACT balloons in 0.5 ml increments.
- If continence is not restored after the addition of 1.0 ml per balloon, the ACT balloons should be imaged to confirm that they are intact and in the correct position.
- Increases in balloon volumes should be stopped once the patient is satisfied with their continence OR both balloons have been filled to their

maximum volume of 8.0 ml each, whichever occurs first.

XIV. Device Explanation

If explanation of one or both devices is required, perform the following:

1. Clean and prepare perineal area.
2. Administer anesthetic as necessary.
3. Palpate device port through skin of labia.
4. Make < 1 cm incision over device port.
5. Expose port through incision.
6. Access port with 0.60 mm (23G) x 28.6 mm needle and syringe. Withdraw all fluid.
7. Applying slow, constant traction to tubing and port, remove device through incision.
8. Close incision with sutures.

XV. Clinical Benefits and Summary of Safety and Clinical Performance (SSCP)

- Patients treated with ACT show a measurable quality of life increase. The majority of patients using the ACT system are either dry or significantly improved.
- The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.
<https://ec.europa.eu/tools/eudamed>

XVI. Device Disposal

- Explanted devices, devices attempted to be implanted, and any devices removed from the sterile packaging should be returned to Uromedica.
- Packaging materials or accessories with a suspected defect should be returned to Uromedica.
- If not returned to Uromedica, explanted devices, devices attempted to be implanted, any devices removed from the sterile packaging, packaging materials, and accessories are to be discarded in accordance with local environmental regulations for biohazard material.

XVII. Clinician Instructions for Completing Patient Implant Card

- Within the shelf box, the ACT implant kit includes a patient information leaflet and an implant card for each individual ACT device included in the kit.
- The implant card(s) are adhered to the inside back cover of the patient information leaflet.
- You are required by EU law (Regulation (EU) 2017/745 on Medical Devices (MDR)) to provide the patient information leaflet and implant card(s) to the patient following implantation.
- Carefully read these clinician instructions prior to filling in the patient's implant card(s).
- Each implant card contains a number of blank fields. Symbols are defined at the end of this document.
- Do not cover or change any of the pre-printed text or symbols on the implant card(s).

NOTE: Repeat the following steps for each implant card included in the patient information leaflet to correspond with the number of ACT devices implanted in the patient. The patient will receive two implant cards for a bi-lateral implant (one for the patient's left side implant and one for the patient's right side implant.)

Step 1: Handwritten information to be filled in by the clinician

- Use a black ballpoint pen. Do not use any other type of pen. Do not use a pencil.
- Use print writing to complete the required information. Do not use cursive writing.
- The clinician must complete the required fields on the FRONT of the implant card.
 - Patient name or patient identification
 - Date of implantation
 - Name and address of the implanting healthcare institution / provider
 - Mark the box indicating the side of the patient the ACT device was implanted on.

Example:

		Indicates patient's left side device
		Indicates patient's right side device

Step 2: Provide the completed implant card(s) and patient information leaflet to the patient.

NOTE: the implant card(s) must not be re-used. If additional ACT implantations are performed, a new implant card(s) must be issued.

Appendix A

Contrast Media	Amount (ml)	+	Sterile H ₂ O (ml)
Cysto-Conray®	8	+	17
Cysto-Conray® II	19	+	6
Iomeron® 200	19	+	6
Iomeron® 250	18	+	7
Iomeron® 300	16	+	9
Iomeron® 350	14	+	11
Iomeron® 400	13	+	12
Isovue® 200	19	+	6
Isovue® 250	15	+	10
Isovue® 300	14	+	11
Isovue® 370	12	+	13
Omnipaque™ 140	23	+	2
Omnipaque™ 180	19	+	6
Omnipaque™ 240	16	+	9
Omnipaque™ 300	13	+	12
Omnipaque™ 350	11	+	14
Ultravist® 150	23	+	2
Ultravist® 240	17	+	8
Ultravist® 300	14	+	11
Ultravist® 370	12	+	13

Cysto-Conray® is a registered trademark of Mallinckrodt Imaging.

Iomeron® and Isovue® are registered trademarks of Bracco.

Ultravist® is a registered trademark of Berlex Laboratories, Inc.

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OMC Medical GmbH