ProACT

Adjustable Continence Therapy for Men
Uromedica, Inc.

Instructions for Use

Please read all information carefully.

All first-time implanting surgeons must have received an in-service training session with a surgeon previously qualified on the ProACT implantation procedure.

Failure to properly follow instructions may result in improper functioning of the device and lead to injury.

Important: This package insert is designed to provide instructions for use of the Uromedica ProACT single-use device and Uromedica reusable surgical instruments. **It is not a comprehensive reference to surgical technique for correcting male urinary incontinence due to intrinsic sphincter deficiency (ISD).** The device should be used only by physicians trained in the surgical treatment of male urinary incontinence due to ISD. These instructions are recommended for general use of the device. Variations in use may occur in specific procedures due to individual technique and patient anatomy.

I. Description of ProACT System

The ProACT therapy consists of:

- A periurethral implantable prosthetic which has a two-lumen extruded tube (2.4 mm outside diameter) with an expandable balloon on the end that is implanted adjacent to the bladder neck. The balloon contains a 1mm radiopaque marker positioned 1mm from the distal tip. Balloon diameter is 25 mm at maximum fill volume (8 ml). The implant contains a subcutaneous titanium port on the end of the tube that is used by the physician for post-operative adjustment.
- A pushwire of 0.8 mm OD and 15 cm in length with a flexible tip on which the ProACT device is mounted and is used to deliver the ProACT through an Implantation Instrument Set.
- A syringe with a non-coring needle.
- Material exposure to user (Implant): 100% silicone.
 - The maximum surface area exposure of this product includes two 9 cm devices, each at maximum fill volume (8 ml). This maximum surface area exposure is 9.015 in².
 - The maximum mass of this product includes two 9 cm devices, each at maximum fill volume (8 ml). This maximum mass is 18.65 g.
- Maximum fill volume of isotonic contrast solution per device: 8 ml.

II. Contents of ProACT International Patient Pack:

Table 1 – ProACT International Model Numbers

Model Number	Device Length	Description
800018-01	12 cm	ProACT Patient Kit: Two (2) prosthetic implants
800018-02	14 cm	ProACT Patient Kit: Two (2) prosthetic implants
800022-01	12 cm	ProACT Revision Kit: One (1) prosthetic implant
800022-02	14 cm	ProACT Revision Kit: One (1) prosthetic implant

All configurations also include:

- One (1) 0.60 mm (23G) x 28.6 mm Huber Needle
- One (1) 5 ml Syringe with Luer Lock Connector
- One (1) Spare Pushwire
- One (1) Instructions for Use Insert
- One (1) Product Label Set

III. Uromedica Surgical Instruments (sold separately)

- A non-sterile Implantation Instrument Set (P/N 750009) consisting of a U-channel sheath and sharp-tipped trocar.
- A non-sterile Blunt Trocar (P/N 750022)

- A non-sterile Tissue Expanding Device (TED II) (P/N 750034)

The Uromedica Blunt Trocar and Uromedica TED are both intended to be used with the U-channel sheath of the Uromedica Implantation Instrument Set. The U-channel sheaths mates with these devices during use, per the instructions below.

IV. Intended Purpose and Indication

Intended Purpose

The ProACT prosthetic implant device is intended, via initial and post-operative volume adjustments, to increase urethral coaptation at the bladder neck or at the apex of the prostate to eliminate unwanted urine loss.

Indication

The ProACT system is designed to provide a permanent implant for the correction of urinary incontinence in adult male patients post prostatectomy and who have either de novo (no prior surgery) or recurrent (failed prior surgery) stress urinary incontinence.

V. Contraindications

- Patients with active systemic or urinary tract infections.
- Patients with incontinence due to detrusor instability.
- Patients with reduced bladder compliance.
- Patients with significant residual volume greater than 100ml after voiding.
- Patients who are presently receiving, or plan to receive within (6) months, radiotherapy and those patients who have received radiotherapy within the last six (6) months.
- Patients suspected of having bladder cancer.
- Patients with unsuccessfully treated bladder stones.
- Patients with hemophilia or bleeding disorders.

VI. Warnings

- The device is an implantable device and is intended for a single use only. Do not re-sterilize. Re-sterilization could result in the mechanical failure of the device and place the patient at undue risk of procedural failure.
- Placement of the ProACT device must be performed by a physician.
- If the physician chooses to use X-ray, adequate attention must be given to the potential radiation exposure.
- Perforation of the bladder, urethra, or rectum during dilation of the tissue from the perineal area to the bladder neck should be avoided.
- Do not use excessive force to advance or withdraw the device if resistance is encountered.
- Do not attempt to place or introduce the device without the supplied pushwire.
- Do not place the device in a location where the patient may have problems with urinary retention or erosion.
- Do not inflate the balloon to more than 8.0 ml.
- Do not place the device in the urethral mucosa. This could result in balloon erosion or migration.
- Use the correct mixture of sterile water and contrast media to provide an isotonic solution for inflation of the ProACT balloons. See attached reference chart (Appendix A).
- Do not use any of the non-sterile accessory surgical instruments without prior cleaning and steam sterilization as indicated in section IX.
- Device wear or puncture may result if the device comes in contact with itself or another device, or from surgical staples, steel suture, or other artifacts from previous interventions present in the region where the device will be positioned. The physician is cautioned to carefully evaluate the implant location.
- The ProACT system is MR-Conditional. A patient with this device can be scanned safely in an MR (Magnetic Resonance Imaging) system immediately after placement under certain conditions as detailed in the

- MRI Information section provided in these Instructions for Use (IFU).
- Device leak may result due to long-term calcification on the surface of the balloon.
- Resolution of adverse events (e.g. device erosion, prosthetic infection, device wear) or suspected device malfunction or changes in its performance may require explant of one or both of the devices. The explanted device(s) can be replaced in 60 to 90 days.

VII. Precautions

- Do not use the ProACT system prior to reading and understanding the Instructions for Use. Product packaging that is damaged or unintentionally opened shall not be used.
- The sterile packaging shall be inspected prior to use.
- Do not attempt to re-sterilize the ProACT device.
- To reduce the risk of needle sticks, always keep the plastic needle cap on the needle until the needle is mated with the syringe and only remove the needle cap prior to using the needle.
 - Proper Disposal: do not attempt to recap a used needle. Do not attempt to bend or break the needle. Immediately dispose of the entire unit (needle, syringe body, and cap) into a puncture-resistant biohazard sharps container.
 - Do not reach into sharps container.
 - Do not leave used needle on the table, even if the cap is fully covering needle.
- Explant devices or devices potentially contaminated should be handled as biohazardous waste.
- The ProACT system should only be used by a physician trained 1) in urologic surgical procedures and 2) who has also completed the in-service training session on ProACT implantation technique conducted by a previously qualified surgeon.
- Patients undergoing a surgical implantation of the ProACT device should be pre-medicated with an antibiotic to reduce the possibility of a postoperative infection.
- Once placed, the ProACT device should be viewed with fluoroscopy visualization, TRUS, or digitally palpated from the anteriolateral rectal wall.
- It is the physicians' responsibility to determine the correct length of the ProACT device.
- Do not implant the ProACT device without the Uromedica, Inc. Implantation Instrument Set and the provided non-coring needle for filling the balloons.
- Do not implant just one ProACT device EXCEPT in the case of bladder or urethral perforation on one side. If a perforation occurs, do not implant a ProACT device on the affected side for at least 60 days. An effective procedure calls for two devices to be placed.
- Improper placement of the ProACT device may lead to a sub-optimal patient clinical response or patient discomfort.
- Maintain the sterile field throughout the procedure.
- If the patient is to undergo any further examination or surgery in the pelvic region after the ProACT device has been implanted, the physician must consider possible adverse effects.
- Do not handle the ProACT device, particularly the balloon, with any sharp instruments or lay it against an abrasive material. Failure of the device may result. Do not grasp the port or strain relief with clamps.
- If the patient is allergic to contrast media, the ProACT balloons should be filled with sterile saline instead of isotonic contrast solution.
- Use the same type of contrast media as used during the original filling of the balloons when performing a post surgical adjustment.
- The Tissue Expanding Device II should not be rotated while the tip is extended. Rotating the tool while the tip is extended may increase the risk of unintended tissue/organ damage.
- The tip of the Tissue Expanding Device II should not be extended towards of the urethra. Extending the tip towards the urethra may increase the risk of perforation.

- The Tissue Expanding Device II should not be used to grasp material. Grasping material may damage the tool.
- Patients with intrinsic sphincter deficiency (ISD) from neurogenic causes may have a worse functional result and may have a higher risk for intra-operative and post-operative complications.
- Patients with prior radiotherapy may have a worse functional result and may have a higher risk of intraoperative perforation, prosthetic migration, and device erosion.

VIII. Potential Adverse Events (Residual Risks):

Potential adverse events from the ProACT system, the attendant equipment and supplies, or the implant procedure include:

Allergic Response
 Anesthetic Risks
 Device Calcification
 Device Malfunction
 Device Wear
 Erosion of the device at the bladder, perineum, rectum, scrotum, or urethra
 Explant of one or both devices
 False Channel Creation
 Hydronephrosis
 Intraoperative Perforation (e.g. bladder, rectum, urethra)
 Pain or Discomfort
 Procedure/adjustment not able to be completed
 Prosthetic Infection
 Prosthetic Migration
 Urinary Complications (e.g. Urinary Difficulty, Frequency, or Urgency, Worsening Incontinence*)
 Urinary Retention
 Urinary Tract Infection
 Vascular Complications (e.g. Bleeding, Hematoma, Induration, Inflammation)
 Worsening Erectile Function
 Wound Dehiscence
 Wound Infection

* "Worsening Incontinence" and "Worsening Erectile Function" includes subjective patient-perceived worsening.
 **Any serious incident associated with the device must be promptly reported to the manufacturer and the relevant competent authority in the Member State where the user or patient is located.

IX. Sterility and Storage

The contents of the undamaged and sealed packaged ProACT device are STERILE. The package has been sterilized with gamma irradiation.

Store product at room temperature with a temperature limit of -30°C to 60°C. Keep the product dry.

The Uromedica Implantation Instrument Set, Blunt Trocar, and Tissue Expanding Device II (TED II) are supplied non-sterile and must be sterilized before initial use and cleaned and sterilized after each use. Follow the reprocessing instructions provided in the Instructions for Use provided with the surgical instruments.

X. MRI Safety Information

MRI Safety Information  MR Conditional	The ProACT, Adjustable Continence Therapy for Men is MR Conditional. A patient with the ProACT, Adjustable Continence Therapy for Men may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.
Parameter	Condition of Use/Information
Static Magnetic Field Strength (T)	1.5-Tesla or 3.0-Tesla
Static Magnetic Field Orientation	Horizontal

Maximum Spatial Field Gradient (T/m and gauss/cm)	50-T/m (5,000-gauss/cm)
RF Excitation Polarization	Circularly Polarized (CP) (i.e., Quadrature-Transmission)
Transmit RF Coil	Any transmit RF coil may be used
Receive-Only RF Coil	Any receive-only RF coil may be used
MR System Operating Mode	Normal Operating Mode
Maximum Whole Body Averaged SAR (W/kg)	2-W/kg (Normal Operation Mode)
Scan Duration and Wait Time	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks)
MR Image Artifact	The presence of this implant produces an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant is located in the area of interest.

XI. Instructions for Implantation

Equipment Requirements

- Image Intensifier and sterile cover
- Transrectal ultrasound and sterile sleeve (Optional for Post Radical Prostatectomy, essential for post Transurethral Resection or enucleation)
- Camera and light source

Instrument Requirements

- Rigid or flexible cystoscopy set up
- Regular Surgical instrument set
- No. 10 blade and handle
- Debakey Forceps
- Rubber shod plain forceps
- Adson Forceps
- Cairns Tissue Forceps
- McIndoe / Metzenbaum Scissors
- Dressing Scissors
- 7" Needleholder
- Uromedica Implantation Instrument Set (sharp trocar and U-channel sheath)
- Uromedica Blunt Trocar
- Uromedica Tissue Expanding Device (TED II)

Accessories

- Jug
- Gallipots x2
- Kidney dish x2
- Bladder syringe
- Syringe
- 3/0 Vicryl rb needle
- 4/0 skin suture
- Foley catheter (typically 14Fr) and drainage bag

Solutions

- One-liter normal sterile saline (0.9%) or sterile water for irrigation during cystoscopy.
- Isotonic contrast solution for filling the device. See Appendix A.
- Sterile saline for soaking the ProACT device prior to implantation.
- Use 50 - 100 mls of contrast and dilute with normal sterile saline (0.9%) to reach a total volume of 250 mls.

Preparing the Patient

- Start oral antibiotic therapy.

- Patient positioned in lithotomy with acute abduction to enable access to perineum. Shave performed at discretion of surgeon.
- Prep perineal and groin region using Betadine. Drape leaving anus exposed if surgeon is utilizing TRUS. Position and drape C arm.
- Perform cystourethroscopy (Solution a). Drain bladder and inflate to 100 - 200ml using bladder syringe and contrast/normal saline mix (Solution d). Cystoscope sheath/obturator retained in place during procedure to visualize and palpate the urethra. Check that bladder neck is identifiable on Image Intensifier (I/I).

Administering Anesthetic

- Administer either general or spinal anesthetic which should be the preference of the urologic surgeon.

Preparing the Contrast Media

- Prepare an isotonic filling solution by diluting a contrast media chosen from Appendix "A" with the prescribed amount of sterile water. The full volumes indicated need not be mixed but the ratios should remain the same so that water does not pass into or out of the device over time due to osmosis. Larger volumes facilitate more accurate mixing.

Prepare the Uromedica Re-usable Surgical Instruments

- The Uromedica Implantation Instrument Set, Blunt Trocar, and Tissue Expanding Device II (TED II) are supplied non-sterile and must be sterilized before initial use, and cleaned and sterilized after each use. Follow the reprocessing instructions provided in the Instructions for Use provided with the surgical instruments.

Implantation of the ProACT device in Post Radical Prostatectomy Patients

- Two transverse perineal incisions are made about 1 cm long and 1 cm posterior to base of the scrotum, immediately posterior to the urethra on either side of the midline.
- Metzenbaum / McIndoe scissors may be used to dissect adipose and muscle layers on each side of the midline.
- Using finger palpation (or Metzenbaum Scissors) through the incision, identify the tip of the fascial triangle created between the urethra and inferior pubic rami. Using a pair of artery forceps, open the incision large enough to allow the introduction of the Implantation Instrument Set (sharp trocar mated to U-channel sheath). The Implantation Instrument Set should be directed toward the inferior pubic rami and should push through the pelvic floor.
- Advance the Implantation Instrument Set and, using controlled force and a twisting motion (rather than pushing), dilate proximally towards the bladder neck. Retain the cystoscope in a horizontal position and use it as a guide by directing the Implantation Instrument Set parallel in a superolateral plane to the membranous urethra. Use fluoroscopy frequently to observe position and avoid bladder or urethral perforation or being more lateral to the urethra than 5 mm. Given the expected presence of scar tissue around the bladder neck, manipulate the Implantation Instrument Set slowly, using only a twisting motion to prevent bladder perforation. The sharp trocar may be replaced by the blunt trocar (mated to the U-channel sheath) to advance the last 0.5 to 1.0 cm up to the bladder neck to reduce the likelihood of perforation of the bladder or urethra. Extra care should be taken in reaching the last 0.5 cm to 1 cm distal to the bladder neck, as this is where bladder perforation is most likely to occur. Once the tip of the trocar appears to be at the bladder neck, stop and remove the trocar from the U-channel sheath.

NOTE: Urethral or bladder perforation can be identified when the trocar is removed and urine is lost out of the

- U-channel sheath. In the event that a urethral perforation occurs, the U-channel sheath should be withdrawn gently, and a balloon should not be implanted on the perforated side. Once recognized, the perforation can be treated by the placement of a urethral Foley catheter for three (3) days.
5. On removal of the trocar, it is important to gently advance the U-Channel sheath proximally to occupy the space created by the tip of the trocar. At this point, if the patient has fibrotic tissue, the Tissue Expanding Device (TED) may be used to dilate the tissue, creating a space for the ProACT balloon. The TED can be locked within the U-Channel sheath. The TED handles are then squeezed together to open the jaw, facilitating tissue dilation. The TED handles must be opened fully (which closes the jaw fully) prior to removing the TED from the U-Channel sheath.
 6. Choose a ProACT device with the appropriate length. A 14 cm length ProACT device is chosen most often. Centimeter markings on the trocar may aid in choosing the proper length ProACT device. Once the tip of the trocar is at the incision site, count the number of centimeters marked on the back side of the trocar and add 4 cm for the final approximate length of the ProACT device to accommodate the positioning of the ProACT port from the incision.
 7. Prime the devices. Remove device from tray and draw the prepared isotonic filling solution (Solution B) into the syringe provided. Attach needle (included) to syringe and inflate the device up to 1.5 ml via the port. Hold the syringe with the plunger towards the ceiling and withdraw solution and any air bubbles. Micro bubbles will not affect device. Aim to remove all the air, leaving a fluid lining within the ProACT device. It should form a tri-wing configuration to facilitate insertion. It may be necessary to form the balloon into a tri-wing configuration with your fingers as the fluid is aspirated. Once the device is primed, remove the needle from the port and submerge device in sterile saline (Solution C).
 8. Remove one device from the sterile saline and ensure that the pushwire is pushed to the tip of the balloon. Coat the balloon with surgical lubricant and place one "wing" directly into the channel of the U-channel sheath.
 9. Advance the ProACT device through the U-channel sheath to the selected position using the pushwire. Check that the balloon is at the bladder neck using fluoroscopy.
 10. Once the balloon is in the correct position, pull back the U-Channel sheath approximately 2 cm to allow room for the balloon to inflate. Prior to inflating the balloon, check to see that the balloon has not moved. Inject up to 1.5 ml of Solution B into the device port using the provided needle. Check the position and inflation using fluoroscopy again.
 11. Leave the pushwire in place. Repeat the previous steps on the contralateral side unless performing a unilateral revision procedure.
 12. Once the second balloon is in place, use fluoroscopy to check that both balloons are positioned correctly. If one or both balloons are too distal to the bladder neck, deflate the balloon, remove the device (soak in sterile saline), and redilate the same tract using the Implantation Instrument Set to place at a more proximal position to the bladder neck. If the balloons are in the correct position, remove the push-wires.
 13. Confirm the absence of urethral perforation by cystoscopic examination. If a urethral perforation has occurred, the balloon on the side of the perforation must be removed (see precautions and warnings).
 14. Using artery forceps, create a pocket passing from the initial incision into the postero-lateral wall of the scrotum to place the ports. Ports should be placed under the skin in the scrotum so that they will not cause discomfort to the patient when seated and should be superficial enough to permit easy access for adjustments.
NOTE: Rubber shod clamps or forceps should be used when grasping the ports to avoid tearing the silicone.
 15. Close the incision in two layers using resorbable suture for the subcutaneous tissue and the skin.
 16. The patient should continue to take the same oral antibiotic given preoperatively.

17. Confirm the patient can void spontaneously. If the patient cannot void spontaneously, place a Foley catheter until anaesthesia effects have subsided, then have the patient void. If the patient is still unable to void, decrease the balloon volume in 0.5 ml increments from both balloons until the patient can spontaneously void.

Implantation of the ProACT device in patients who have undergone Transurethral resection or enucleation for Benign Prostatic Hyperplasia

1. Two transverse perineal incisions are made about 1 cm long and 1 cm posterior to base of the scrotum, immediately posterior to the urethra on either side of the midline.
2. Metzenbaum / McIndoe scissors may be used to dissect adipose and muscle layers on each side of the midline.
3. Using finger palpation (or Metzenbaum Scissors) through the incision, identify the tip of the fascial triangle created between the urethra and inferior pubic rami. Using a pair of artery forceps, open the incision large enough to allow the introduction of the Implantation Instrument Set (sharp trocar mated to U-channel sheath). The Implantation Instrument Set should be directed toward the inferior pubic rami and should push through the pelvic floor.
4. Advance the Implantation Instrument Set and, using controlled force and a twisting motion (rather than pushing), dilate proximally towards the bladder neck. Retain the cystoscope in a horizontal position and use it as a guide by directing the Implantation Instrument Set parallel in a superolateral plane to the membranous urethra. Care should be taken to avoid puncturing the prostate. Use fluoroscopy frequently to observe position and avoid bladder or urethral perforation or being more lateral than 5mm or too medially into the urethra. The absence of scar tissue should allow for the creation of a space at the level of the membranous urethra, lateral to the external sphincter and dorsolateral to the apex of the prostatic remnant. The use of fluoroscopy can assist in determining the length of the prostatic urethra in order to estimate anatomic position of the apex more easily. The measurements should be taken when the trocar is stationary, as the prostate and bladder base will move when prodded by the trocar. Once the tip of the sharp trocar appears to be at the level of the apex, stop and remove the sharp trocar. Extra care should be taken in reaching the last 0.5 cm to 1 cm distal to the apex of the prostate, as this is where perforation is most likely to occur. The blunt trocar can be used to advance the last 0.5 to 1.0 cm up to the apex of the prostate to reduce the likelihood of perforation.
NOTE: Bladder or urethral perforation can be easily identified when the trocar is removed and urine is lost out of the U-channel sheath. In the event that this occurs, the U-channel sheath should be withdrawn gently, and the ProACT device should not be implanted until perforation is healed. Once recognized, the bladder or urethral perforation can be treated by the placement of a urethral Foley catheter for three (3) days.
5. On removal of the trocar, it is important to gently advance the U-Channel sheath proximally to occupy the space created by the tip of the trocar. At this point, if the patient has fibrotic tissue, the Tissue Expanding Device (TED) may be used to dilate the tissue, creating a space for the ProACT balloon. The TED can be locked within the U-Channel sheath. The TED handles are then squeezed together to open the jaw, facilitating tissue dilation. The TED handles must be opened fully (which closes the jaw fully) prior to removing the TED from the U-Channel sheath.
6. Choose a ProACT device with the appropriate length. A 14 cm length ProACT device is chosen most often. Centimeter markings on the trocar may aid in choosing the proper length ProACT device. Once the tip of the trocar is at the incision site, count the number of centimeters marked on the back side of the trocar and add 4 cm for the final approximate length of the

ProACT device to accommodate the positioning of the ProACT port from the incision.

7. Prime the devices. Remove device from tray and draw the prepared isotonic filling solution (Solution B) into the syringe provided. Attach needle (included) to syringe and inflate the device up to 1.5 ml via the port. Hold the syringe with the plunger towards the ceiling and withdraw solution and any air bubbles. Micro bubbles will not affect device. Aim to remove all the air, leaving a fluid lining within the ProACT device. It should form a tri-wing configuration to facilitate insertion. It may be necessary to form the balloon into a tri-wing configuration with your fingers as the fluid is aspirated. Once the device is primed, remove the needle from the port and submerge device in sterile saline (Solution C).
8. Remove one device from the sterile saline and ensure that the pushwire is pushed to the tip of the balloon. Coat the balloon with surgical lubricant and place one 'wing' directly into the channel of the U-channel sheath.
9. Advance the ProACT device through the U-channel sheath to the selected position using the pushwire. Check that the balloon is at the correct position using fluoroscopy.
10. Once the balloon is in the correct position, pull back the U-Channel sheath approximately 2 cm to allow room for the balloon to inflate. Prior to inflating the balloon, check to see that the balloon has not moved. Inject up to 1.5 ml of Solution B into the device port using the provided needle. Check the position and inflation using fluoroscopy again.
11. Leave the pushwire in place. Repeat the previous steps on the contralateral side unless performing a unilateral revision procedure.
12. Once the second balloon is in place, use fluoroscopy to check that both balloons are positioned correctly. If one or both balloons are judged to be too distal to the apex of the prostatic remnant, deflate the balloon, remove the device (soak in sterile saline), and redilate the same tract using the Implantation Instrument Set to place at a more proximal position to the apex of the prostatic remnant. If the balloons are in the correct position, remove the push-wires.
13. Confirm the absence of urethral perforation by cystoscopic examination. If a urethral perforation has occurred, the balloon on the side of the perforation must be removed (see precautions and warnings).
14. Using artery forceps, create a pocket passing from the initial incision into the postero-lateral wall of the scrotum to place the ports. Ports should be placed under the skin in the scrotum so that they will not cause discomfort to the patient when seated and should be superficial enough to permit easy access for adjustments.
NOTE: Rubber shod clamps or forceps should be used when grasping the ports to avoid tearing the silicone.
15. Close the incision in two layers using resorbable suture for the subcutaneous tissue the skin.
16. The patient should continue to take the same oral antibiotic given preoperatively.
17. Confirm the patient can void spontaneously. If the patient cannot void spontaneously, place a Foley catheter until anaesthesia effects have subsided, then have the patient void. If the patient is still unable to void, decrease the balloon volume in 0.5 ml increments from both balloons until the patient can spontaneously void.

XII. Post-Operative Care

- The patient is usually discharged on the day of the surgery. Ensure that the patient can void spontaneously prior to discharge from the surgical facility. Most surgeons choose to send their patients home with a course of antibiotics. The patient should be advised on the use of absorbent pads or condom catheters until the desired level of continence is achieved. Adjustments may begin six weeks after the ProACT implantation procedure. Inform the patient that multiple adjustments may be required to achieve the desired level of continence.
- Provide patient education about the operation. The patient should be advised to avoid undue

strenuous activity immediately post-operatively as directed by the physician. Failure to do so could lead to device migration. Instruct the patient to contact the physician if they experience redness, swelling, or drainage of the incision site, experience a return of incontinence symptoms, or are unable to urinate. Provide the patient with information on the MRI conditional status of the ProACT implant and the impact that ProACT might have on travel security screenings (titanium and radiopaque components).

- After the postoperative healing period, the surgeon should maintain contact with the patient to evaluate the function of the device. During an evaluation, the physician should ask the patient about how ProACT is functioning and if the patient has noticed any changes in the function. If the patient's incontinence symptoms have returned, or there is infection or erosion present, revision or removal surgery may be necessary. In the event of a revision surgery, follow the same explanation and implantation techniques outlined elsewhere in this document.

XIII. Post-Operative Adjustment of Patient

- After edema from the surgery has subsided (approximately six weeks postoperatively) the patient may become incontinent, in which case, adjustments should be performed.
- Adjustments may be conducted in the office. Topical anesthetic may be used, if needed to minimize patient discomfort during palpation of the port and insertion of the needle.
- Do not increase the balloon volume until six weeks post-implantation. Increases in balloon volume performed before six weeks post-implantation are not advised and can result in device migration.
- Use a 0.60 mm (23G) x 28.6 mm needle for post operative balloon volume adjustments. Use of other needles may result in damage to the device port.
- Clean and prepare skin over the adjustment port. Palpate the implanted port to allow percutaneous adjustment of balloon volume. Adjust the balloon volume by accessing the port using a 0.60 mm (23G) x 28.6 mm needle. The ProACT balloon volume should be increased in 0.5 ml increments. To maintain isotonicity, use the same inflation media (diluted contrast media or sterile saline) when performing a balloon volume adjustment as was used during the initial inflation. CAUTION: If resistance is felt, do not continue pushing the needle into the port. Pull the needle back and ensure device tubing is not kinked. Pushing the needle into the edge of the port can cause damage to the port, which may require explant of the device.
- Make sure that patient can spontaneously void prior to leaving the office.
- If the patient experiences urinary retention, remove volume from each of the ProACT balloons in 0.5 ml increments.
- If continence is not restored after the addition of 1.0 ml per balloon, the ProACT balloons should be imaged to confirm that they are intact and in the correct position.
- Increases in balloon volumes should be stopped once the patient is satisfied with their continence OR both balloons have been filled to their maximum volume of 8.0 ml each, whichever occurs first.

XIV. Device Explanation

If explantation of one or both devices is required, perform the following:

1. Clean and prepare perineal area.
2. Administer anesthetic as necessary.
3. Palpate device tubing through skin of perineum.
4. Make < 1 cm incision over device tubing.
5. Expose port through incision.

6. Access port with 0.60 mm (23G) x 28.6 mm needle and syringe. Withdraw all fluid.
7. Applying slow, constant traction to tubing and port, remove device through incision.
8. Close incision with sutures.

XV. Clinical Benefits and Summary of Safety and Clinical Performance (SSCP)

- Patients treated with ProACT show a measurable quality of life increase. The majority of patients using the ProACT are either dry or significantly improved.
- The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

<https://ec.europa.eu/tools/eudamed>

XVI. Device Disposal

- Explanted devices, devices attempted to be implanted, and any devices removed from the sterile packaging should be returned to Uromedica.
- Packaging materials or accessories with a suspected defect should be returned to Uromedica.
- If not returned to Uromedica, explanted devices, devices attempted to be implanted, any devices removed from the sterile packaging, packaging materials, and accessories are to be discarded in accordance with local environmental regulations for biohazard material.

XVII. Clinician Instructions for Completing Patient Implant Card

- Within the shelf box, the ProACT implant kit includes a patient information leaflet and an implant card for each individual ProACT device included in the kit.
- The implant card(s) are adhered to the inside back cover of the patient information leaflet.
- You are required by EU law (Regulation (EU) 2017/745 on Medical Devices (MDR)) to provide the patient information leaflet and implant card(s) to the patient following implantation.
- Carefully read these clinician instructions prior to filling in the patient's implant card(s).
- Each implant card contains a number of blank fields. Symbols are defined at the end of this document.
- Do not cover or change any of the pre-printed text or symbols on the implant card(s).

NOTE: Repeat the following steps for each implant card included in the patient information leaflet to correspond with the number of ProACT devices implanted in the patient. The patient will receive two implant cards for a bi-lateral implant (one for the patient's left side implant and one for the patient's right side implant.)

Step 1: Handwritten information to be filled in by the clinician

- Use a black ballpoint pen. Do not use any other type of pen. Do not use a pencil.
- Use print writing to complete the required information. Do not use cursive writing.
- The clinician must complete the required fields on the FRONT of the implant card.
 - Patient name or patient identification
 - Date of implantation
 - Name and address of the implanting healthcare institution / provider
 - Mark the box indicating the side of the patient the ProACT device was implanted on.

Example:

	Indicates patient's left side device
	Indicates patient's right side device

Step 2: Provide the completed implant card(s) and patient information leaflet to the patient.

NOTE: the implant card(s) must not be re-used. If additional ProACT implantations are performed, a new implant card(s) must be issued.

Appendix A

Contrast Media	Amount	+	Sterile H ₂ O
Cysto-Conray®	8ml	+	17ml
Cysto-Conray® II	19ml	+	6ml
Iomeron® 200	19ml	+	6ml
Iomeron® 250	18ml	+	7ml
Iomeron® 300	16ml	+	9ml
Iomeron® 350	14ml	+	11ml
Iomeron® 400	13ml	+	12ml
Isovue® 200	19ml	+	6ml
Isovue® 250	15ml	+	10ml
Isovue® 300	14ml	+	11ml
Isovue® 370	12ml	+	13ml
Omnipaque™ 140	23ml	+	2ml
Omnipaque™ 180	19ml	+	6ml
Omnipaque™ 240	16ml	+	9ml
Omnipaque™ 300	13ml	+	12ml
Omnipaque™ 350	11ml	+	14ml
Ultravist® 150	23ml	+	2ml
Ultravist® 240	17ml	+	8ml
Ultravist® 300	14ml	+	11ml
Ultravist® 370	12ml	+	13ml

Cysto-Conray® is a registered trademark of Mallinckrodt Imaging.
Iomeron® and Isovue® are registered trademarks of Bracco.
Ultravist® is a registered trademark of Berlex Laboratories, Inc.

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OMC Medical GmbH